

Sakhi — Anaemia Screen Test Card

Hold ONE colour block up to the webcam, filling the pink circle, in good light.
It mimics eyelid tissue so you can demo the screen without using your eye.

A. Healthy conjunctiva → expect: low pallor (none/mild)



B. Pale / possible anaemia → expect: higher pallor (moderate/strong)



Tips:

- Bright, even light. Avoid glare and coloured lamps.
- Fill the circle with the block; hold steady until the light turns green/amber.
- This is a rough demo aid, not a medical calibration.
- MediaPipe won't see an 'eye' in a flat card, so if the eye-gate blocks capture, test on a phone where the gate falls back, or point at a real lower eyelid.